



TIKI SIZE GUIDE EN 1091-61

Size guide

Size guide for the Tiki medical. For guidance only.

Measure the distance between your upper forehead and the lower edge of your chin.

< 190 mm	Small
190–200 mm	Medium
> 200 mm	Large

Important! Before using the Tiki respirator you have to do a seal check according to the instructions in the manual.



TIKI SAFETY AB • SKÅLLSTAVÄGEN 9 • SE-19740
BRO • SWEDEN • TIKIMED.COM
MADE IN SWEDEN



The complete mask is certified by NB according to ISO EN 12942.



The exhalation filter is a Medical device Class 1 product intended to be used in combination with the complete mask.